



Amber ALE



GLUTEN FREE RECIPE KIT

CUT THIS OUT FOR
YOUR TAP HANDLE

AMBER ALE

Gluten Free Recipe Kit



EXTRACT (KIT910)

6 lbs Sorghum Extract

1 lb Candi Syrup D180



HOPS

1 oz Kent Goldings – Boil 60 min

0.5 oz Northern Brewer – Boil Last 30 min



OTHER ADDITIONS

Clarifier – Boil last 5 min

Dry Yeast Packet – Rehydrate & pitch into fermentation vessel



FITS OUR TAP HANDLE D1282 PERFECTLY



AMBER ALE

Gluten Free Recipe Kit

RECIPE TIPS

To hydrate dry yeast, sprinkle yeast in 4 oz of distilled water or wort at 77–84°F. Do not stir. Let rest for 15 to 30 minutes. Begin stirring occasionally for an additional 30 minutes. Pitch the resultant cream into your fermentation vessel.

We recommend a fermentation temperature in the 66–68°F range.

NEED HELP?

SEARCH

AMBER ALE
GLUTEN FREE KIT
on MoreBeer.com
to see detailed brewing
instructions & descriptions

LEARN MORE

Over 250 Articles at
morebeer.com/articles
&
Over 125 Videos at
morebeer.com/videos

BREW DAY QUESTIONS?

Our friendly staff
are ready to help!

GIVE US A CALL
800-600-0033
OR EMAIL US 24/7
info@morebeer.com

BEER STATS

EST. ORIGINAL GRAVITY
1.045-51

ACTUAL O.G.

SRM: 16.7

IBU'S: 23.7

EST. ABV %: 5.3

ACTUAL ABV %:

DATE BREWED:

NOTES:

NOTES



MY RATING:

1 2 3 4 5

TASTING DATE: _____

APPEARANCE

AROMA

TASTE

OVERALL

REVIEW THIS KIT ON MOREBEER.COM

RECORD KEEPING



DATE BREWED _____

GALLONS IN FERMENTER _____

ORIGINAL GRAVITY _____

TEMP OF WORT AT PITCH _____

LAG TIME _____

FERMENTATION TEMP _____

DAYS IN FERMENTER _____

FINAL GRAVITY _____

YEAST STRAIN _____

FOR ALL GRAIN BREWERS



STRIKE WATER TEMP _____

MASH TEMP _____

BREWING CALCULATORS

morebeer.com/calc

Use our free online Brewing Calculators to calculate Alcohol %, Boil Off %, Hydrometer Correction, & More!